

# Seg-WAY<sup>TM</sup> EPF

## Endoscopic Plantar Fasciitis Release (EPF)

Our Unique Seg-WAY Endoscopic System gives you greater control and unparalleled visualization. It delivers the tactile feel of an endoscopic scalpel in a minimally invasive single-portal medial approach that can help reduce the lengthy recovery time associated with these types of surgeries. Unlike most competitive systems, the scope functions independently from the knife, rasp, and probe providing ultimate system flexibility.

# Seg-WAY™ EPF



Independent retrograde knife operation



Left and Right guide for blade handling preference




Probe for safe release



Rasp for better visualization before release

## Features and Benefits

- Uni-Minimal Incision - 8mm–10mm
- Minimal Tissue Damage - Insertion of guide as far as needed
- Faster Surgery Time - Approximately 2–5 minutes
- Potential for patient to be mobile the next day
- Left and Right Guide - Blade handling preference
- 3x the field of view
- Tactile knife and probe
- System utilizes a standard 4mm 30° scope



Description	Part #
Seg-WAY Instrument Set (Complete)	200-5560
Left Large SS Knife Guide	200-5012
Right Large SS Knife Guide	200-5022
Left Small SS Knife Guide	200-5014
Right Small SS Knife Guide	200-5024
Synovial Dilator / Elevator	200-5040
Rasp, Ligament	200-5050
Probe, Ligament	200-5060
Dilator	000-0327
Single Use Retrograde Knife (62mm)	200-1003
Single Use Retrograde Knife (45mm)	200-1009