

# Seg-WAY<sup>TM</sup>

## ECTR

### Endoscopic Carpal Tunnel Release (ECTR)

Seg-WAY is the first ECTR system designed to properly position the incision in the Ulnar Safe Zone. It is also the only system that is anatomy and patient specific. It provides a level of safety and quality that is unavailable using any other method. The system provides the ability to rasp and probe prior to release for enhanced visualization and safety.

# Seg-WAY<sup>TM</sup> ECTR

## Features and Benefits

### SAFETY FIRST

- Anatomy Specific - Left and right guides place the blade in the Ulnar Safe Zone
- Larger Field of View - a wider, longer, panoramic field of view
- Rasp - to visualize transverse fibers before release

### VERSATILITY

- Patient Specific - Different size guides for smaller and larger hands

### ASSURANCE

- Visually probe the distal end before release

### CONFIDENCE

- Probe fibers to ensure complete release

### UNIPORTAL

- Smaller incision in the wrist



Anatomy and patient specific left and right guide system




Independent retrograde knife operation



Probe distal end and uncut fibers to ensure a complete release



Rasp to visualize transverse fibers before release



Description	Part #
Seg-WAY Instrument Set (Complete)	200-5560
Left Large SS Knife Guide	200-5012
Right Large SS Knife Guide	200-5022
Left Small SS Knife Guide	200-5014
Right Small SS Knife Guide	200-5024
Synovial Dilator / Elevator	200-5040
Rasp, Ligament	200-5050
Probe, Ligament	200-5060
Dilator	000-0327
Retrograde Ligament Knife (Single Use)	200-1003



40 General Warren Blvd, Suite 100 | Malvern, PA 19355 U.S.A.  
610.989.8080 | [www.tricemedical.com](http://www.tricemedical.com) | [customerservice@tricemedical.com](mailto:customerservice@tricemedical.com)