


# Seg-WAY<sup>TM</sup> ECTR-d



## Disposable Endoscopic Carpal Tunnel Release Kit

Seg-WAY ECTR-d is the first, fully disposable endoscopic carpal tunnel release system on the market designed to properly position the incision in the Ulnar Safe Zone. The new, innovative system provides enhanced visualization and safety while enabling the procedure to be performed under local anesthetic. The tray and instruments are completely disposable, resulting in quicker setup and turnaround times.



mi-eye<sup>TM</sup>

# Seg-WAY<sup>TM</sup> ECTR-d

## Features & Benefits:

### EFFICIENCY

- Entire tray and all instruments are fully disposable
- Quick turnaround time; more procedures
- Can be done in-office with local anesthetic

### SAFETY FIRST

- Anatomy Specific - Left and right guides place the blade in the Ulnar Safe Zone
- Rasp - to visualize transverse fibers before release
- Probe - to visualize and locate ligament distal end before release and locate uncut fibers to ensure complete release

### VERSATILITY

- Compatible with 2.3 mm or smaller diameter off-axis angled endoscopes



Tenotomy Scissors




Right/Left Guides



Adson Tissue Forceps



Scalpel



Retractor



Retrograde Knife



Dilator/Elevator



Probe/Rasp



Description	Part #
Disposable Right Guide	1-10-0462
Disposable Left Guide	1-10-0463
Dilator/Elevator	1-10-0470
Disposable Probe/Rasp	1-10-0456
Retrograde Ligament Knife	1-10-0430
Adson Tissue Forceps	2-10-0256
Tenotomy Scissors	2-10-0257
Disposable Scalpel	2-10-0258
Ragnell Retractor	2-10-0262