

Seg-WAYTM EGR

Endoscopic Gastrocnemius Recession (EGR)

Our Unique Seg-WAY Endoscopic System gives you greater control and unparalleled visualization. It delivers the tactile feel of an endoscopic scalpel in a minimally invasive approach that can help reduce the lengthy recovery time associated with these types of surgeries. Unlike most competitive systems, the scope functions independently from the knife, rasp, and probe providing ultimate system flexibility.

Seg-WAY™ EGR



Independent retrograde
knife operation



Left and Right guide for blade handling
preference and measurement




Probe distal end and uncut fibers
to ensure a complete release



Rasp to visualize transverse
fibers before release

Features and Benefits

- Uni-Minimal Incision - 8mm–10mm
- Minimal Tissue Damage
- Faster Surgery Time - Approximately 2–5 minutes
- Patient up and walking the next day
- Left and Right Guide - Blade handling preference
- 3x the field of view
- Tactile knife and probe
- System utilizes a standard 4mm 30° scope



Description	Part #
Seg-WAY Instrument Set (Complete)	200-5560
Left Large SS Knife Guide	200-5012
Right Large SS Knife Guide	200-5022
Left Small SS Knife Guide	200-5014
Right Small SS Knife Guide	200-5024
Synovial Dilator / Elevator	200-5040
Rasp, Ligament	200-5050
Probe, Ligament	200-5060
Dilator	000-0327
Retrograde Ligament Knife (Single Use)	200-1003



40 General Warren Blvd, Suite 100 | Malvern, PA 19355 U.S.A.
610.989.8080 | www.tricemedical.com | customerservice@tricemedical.com